

Any clinician with experience treating hospitalized patients deals with hypoxic respiratory failure on a daily basis. For the respiratory specialist there is a frequent challenge that occurs when these patients on supplemental oxygen by facemask also require nebulized inhaled therapy. The SINODA device is a simple and elegant solution to this all too common problem that allows for uninterrupted flow of oxygen while allowing simultaneous inhaled medication administration. The fact that it is compatible with multiple types of oxygen interfaces, including tracheostomy collars, makes the SINODA device seemingly indispensable.

A handwritten signature in black ink, consisting of a series of loops and a long horizontal stroke extending to the right.

Andrew Toma, MD

Board Certified Pulmonary Medicine and Critical Care Medicine