


SINODA™ device set

I was given an opportunity to evaluate and review the SINODA™ device set. For patients who need combination therapy, SINODA™ is a flexible, multipurpose tool. It is easy to use and works with most systems, so specific adapters are not necessary when treatment speed is crucial. Clinicians can perform diagnostic MRIs without removing vital breathing support because there are no metal or magnetic components. Reduces cross-contamination and nosocomial ICU infections. Easy to use, restock, and dispose of. Physical slots automatically expand and contract in response to variations in body and room temperature to prevent tight lockage.

The three-piece modular respiratory adapter known as the SINODA™ device set enables caregivers to administer continuous oxygen or humidification therapy in addition to aerosolized medication. To facilitate various connection configurations, this one-time connector is divided into three segments that can be oriented in any direction. This connector works with face masks, T-collars, and current gas and aerosolized medicine delivery systems.

Thank you,

A handwritten signature in black ink, appearing to read 'Irene Valdez-Badere', with a long horizontal flourish extending to the right.

Irene Valdez-Badere, RRT, RN, MSN